



Axygen® Mini Plate Spinner

- ▶ Use before and after thermal cycling to improve PCR yield
- ▶ Quickly spins down droplets and condensation
- ▶ Maximum capacity: 2 standard PCR or microtiter microplates
- ▶ Accepts skirted and non-skirted and all standard PCR plates
- ▶ Ambient operating range 4°C to 35°C
- ▶ Exceptionally small footprint (16 x 16 x 18 cm)
- ▶ Starts and stops in seconds



Specifications

Speed	2,500 rpm
Maximum Capacity	2 standard PCR microplates
Ambient Operating Range	4°C to 35°C
Dimensions	6 x 6 x 7 in. (16 x 16 x 18 cm)



Ordering Information

Cat. No.	Description
PLATESPINNER-230EU	Axygen mini plate spinner centrifuge, 230V-EU plug

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences, or contact our Scientific Support Team at ScientificSupportEMEA@corning.com.

Corning Incorporated Life Sciences Europe

Corning BV
Fogostraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 (0) 20 659 60 51
Fax: +31 (0) 20 659 76 73
cceurnl@corning.com
www.corning.com/lifesciences

Support Offices

Eu ROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

Distributed by :

Comercial Hospitalaria Grupo 3

S.L. Rúa do Arroncal, nº9, Vial C,
Nave 4 C Parque Empresarial
Porto do Molle
36350 Nigrán
www.chgrupo3.com

CORNING | **FALCON**® **AXYGEN** **GOSSELIN**™

Offer valid in Spain only and void where prohibited by government, law or company policy. Offer valid for end-user customers only for purchases until December 31, 2017 or while supply lasts.

For a listing of trademarks, visit www.corning.com/clstrademarks.
All other trademarks are the property of their respective owners.