

BeneHeart D1 Quick Reference Guide
BeneHeart D1快速操作卡

 **IF PATIENT IS UNRESPONSIVE
AND NOT BREATHING NORMALLY**
如果病人无反应且无正常呼吸


1

PRESS "ON"
BUTTON
按下开机按键



2

APPLY PADS
FOLLOW AED
INSTRUCTIONS
按照AED指示
粘贴电极片



3

IF INSTRUCTED
PRESS "SHOCK"
BUTTON
如有指示
按下电击按键

